




MiDICOM



Medical Imaging Systems



DICOM Converter Solution

Microtek's MiDICOM provides a reliable converter solution for bringing texts and images together to create DICOM-format medical reports. By applying with MiDICOM, various e-files formats, such as Microsoft Word, Excel, Powerpoint, Adobe PDF, and JPEG, can be transformed into DICOM format in seconds. These DICOM-converted files can be directly previewed and uploaded to PACS/RIS systems, quickly and efficiently managed and shared among the medical professionals.

In addition, working MiDICOM together with Microtek's scanners, users can tremendously reduce paper and chemical waste associated with standard physical modalities, which means no more fuzzies to deal with physicians' handwritings, diagnostic reports or any hard-copy files. All files will be electronicized into DICOM format quickly.

Features

- Downloads data directly from the worklist server
- Allows to create patient information manually
- Converts various file formats easily into DICOM format
- Uploads data to PACS/RIS server

Convert and Upload Easily to PACS/RIS

1. Select "MiDICOM"
2. Print to covert file
3. Send to PACS/RIS


Download patient data




Worklist
Server



Non DICOM format



DICOM 3.0



MiDICOM

Scan Directly to DICOM Format

Working MiDICOM together with Microtek's scanners, users can easily convert any paper document to DICOM format files.



FS-1600XL




MiDICOM



ASDI-4020



Microtek International Inc.
No. 6, Industry East Road 3,
Science-Based Industrial Park,
Hsinchu, Taiwan 30077



EVESAR GmbH,
Linsellesstraße 127-129
D-47877 Willich, Germany

MICROTEK
DICOM READY

www.microtek.com