



MiPAX

Web-Based PACS

Medical Imaging Systems



Features

- CE MDD certified
- Supports telemedicine by Web-based PACS
- Direct access through internet to review medical records, preview images, and download DICOM files
- Supports JPEG2000 lossless and lossy compression format
- On-line review medical images, and instant messages through E-Board function
- Easy to use and effective to communicate
- Offers better image quality, performance, and secure e-medical records



User Interface of MiPAX-WebView

Medical Specialists can direct access through internet to review medical records, preview image, and download the files in JPEG2000 lossless and lossy compression format.



Direct Access through E-Board

On-line Communication for Review and Notations from Multiple Sites



E-Board

Medical Specialists of multiple sites can easily make on-line conference meeting, review medical images, and input messages instantly through E-Board of MiPAX-WebView.



Specifications of Optional X-Ray Digitizer

Product Name	Medi-6000 Plus
Model Type	MSS-600XU2L
Scanning Area	2.5" x 2.5" (Min.) – 14" x 35"(Max.)
Scanning Speed	12 seconds @ 300 dpi, 14" x 17"
Scanning Mode	8 / 12 / 16-bit grayscale
Optical Resolution	600 dpi , (42 µm spot size, 12 lp/mm)
Dmax	4.0
Interface	Hi-Speed USB (USB 2.0)
Dimensions (L x W x H)	305.64 x 474 x 234.5 mm (w/o tray)
Weight	8.7 kg
Power Supply	100 to 240 VAC, 47 ~ 63 Hz
Power Consumption	54.9 W (Max.)
Certifications	FDA, SFDA, CCC, CE, FCC, RoHS
System Requirements	- Microsoft Windows® XP Professional SP3 - CPU : Intel® Core2 Duo E7400 - RAM : 2GB DDR2 800 - Hard Disk : 250GB - Gigabit Ethernet LAN - Monitor : 19" (4:3 / 1280x1024) - DVD-RW



Product Name	Medi-7000
Model Type	MSS-600XU2L
Scanning Area	2.5" x 2.5" (Min.) – 14" x 35"(Max.)
Scanning Speed	7 seconds @ 300 dpi, 14" x 17"
Scanning Mode	8 / 12 / 16-bit grayscale
Optical Resolution	600 dpi , (42 µm spot size, 12 lp/mm)
Dmax	4.0
Interface	Hi-Speed USB (USB 2.0)
Dimensions (L x W x H)	330 x 526 x 435 mm (w/o tray) 738 x 526 x 696 mm (w/tray)
Weight	17.5 kg
Power Supply	100 to 240 VAC, 47 ~ 63 Hz
Power Consumption	54.9 W (Max.)
Certifications	FDA, SFDA, CCC, CE, FCC, RoHS
System Requirements	- Microsoft Windows® XP Professional SP3 - CPU : Intel® Core2 Duo E7400 - RAM : 2GB DDR2 800 - Hard Disk : 250GB - Gigabit Ethernet LAN - Monitor : 19" (4:3 / 1280x1024) - DVD-RW



© 2010 Microtek. Microtek, Medi, MiPAX, ScanWizard are trademarks or registered trademarks of Microtek International, Inc. All other brands, product names and logos herein are trademarks or registered trademarks of their respective holders. Specifications, software bundles and accessories are subject to change without notice. Delivery of technical support services subject to change without notice. Not responsible for typographical errors. Digital Imagery © 2010 Jupiterimages Co. Oct. 31, 2010.

