



# **MiPAX** *Endo* **HD**

**Capture While Recording  
& Link to PACS**

**Full HD**




## Features

- Capable to capture 720p / 1080i / 1080p high quality HD images
- Supports both dynamic video recording and single image capturing
- Capable to connect with Worklist to acquire checking lists of patients, fastening operation time
- Capable to combine captured images together with patients' information, transform them into DICOM format and then upload to PACS directly
- Meets the 3.0 DICOM standard on image capturing, saving and transmitting
- Supports Blue-ray recording in DICOMDIR format and built in CD-Viewer for convenient disk reviewing
- Touch panel brings an easy and intuitive operating way to users



**MICROTEK**  
DICOM READY

## Capture While Recording



## Supported Equipments

Olympus	CV290 CV260 OTV-S190 CV180 OTV-S7Pro
Stryker	1188HD 1288
Pentax	EPK-i7000
Conmed Linvatec	IM4000
Karl Storz	222010 20
da Vinci	da Vinci Si Surgical System

## Powerful Features



For more details on product, please contact your local Microtek sales or dealers.

### Microtek International, Inc.

No. 6 Industry East Road 3, Science-Based Industrial Park, Hsinchu 30075, Taiwan  
Tel: +886-3-577-2155 Fax: +886-3-577-2598 <http://www.microtek.com>

© 2015 Microtek. Microtek, ScanWizard, MiPAX, iPAX, MiDICOM and the MII logo are trademarks or registered trademarks of Microtek International, Inc. All other brands, product names and logos herein are trademarks or registered trademarks of their respective holders. © 2015 123RF Publishing. 20151130

Specifications, software bundles and accessories are subject to change without notice.  
Delivery of technical support services is subject to change without notice.  
Not responsible for typographical errors.

**MICROTEK**  
DICOM READY

[www.microtek.com](http://www.microtek.com)